



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 1 Of 17

PROCESS VALIDATION PROTOCOL

PRODUCT NAME

LENALIDOMIDE CAPSULES 10 MG

STANDARD BATCH SIZE

100000 CAPSULES

PROTOCOL NO.

LD/PVP/004-00



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 2 Of 17

Table of Contents

Sr. No	Content	Page No.
1.	Approval Sheet	3 of 17
2.	Introduction	4 of 17
3.	Objective	4 of 17
4.	Scope	4 of 17
5.	Revalidation	4 of 17
6.	Validation Team	5 of 17
7.	Product and Batch Detail	6 of 17
8.	Equipment & Machine List	7-8 of 17
9.	Batch Weighing Sheet	9 of 17
10.	Approved Raw & Packing List	10 of 17
11.	Process Flow Diagram	11 of 17
12.	Manufacturing Process	12 of 17
13.	Filling & Polishing	13 of 17
14.	Visual Inspection	13 of 17
15.	Sampling Plan	14 of 17
16.	Acceptance Criteria	15 of 17
17.	Deviation	16 of 17
18.	Change Control	16 of 17
19.	Abbreviation	16 of 17
20.	Post Approval	17 of 17



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 3 Of 17

1. Approval Sheet

Protocol Prepared By

Name	Designation	Signature	Date
Sachin Kumar	Sr. Officer QA		03/04/2024

Protocol Checked By

Name	Designation	Signature	Date
Geetesh Kumar	Sr. Executive		03/04/2024
Ajay Bhatia	Sr. Officer QC		03/04/2024

Protocol Approved By

Name	Designation	Signature	Date
N.P. Mishra	Head QA		03/04/2024



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 4 Of 17

2. Introduction:

Three consecutive batches taken for process validation of **Lenalidomide Capsules 10 mg**. If failure in any circumstances further three batches validated till result of all three batches consistent and complying the acceptance criteria.

3. Objective:

The document evidence to provide a high degree of assurance that the manufacturing process of **Lenalidomide Capsules 10 mg** will produce the product meeting its predetermined specification and quality attributes.

4. Scope:

This validation carried out for **Lenalidomide Capsules 10 mg** by monitoring all critical parameters during manufacturing that have impact on product quality.

5. Revalidation:

The manufacturing process is subjected to revalidation when following changes occur.

1. When change in quantity and change in Excipients
2. Change in master formula.
3. Change in manufacturing process.
4. Change in equipment and machinery.
5. Change in vendor of Raw material.



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 5 Of 17

6. Validation team:

Following validation team is responsible for carried out process validation.

Sr. No.	Designation	Department	Responsibility
1.	IPQA	Quality Assurance	Sampling of sample as per sampling plan
2.	Maintenance Head	Maintenance	For smooth working of utilities
3.	Sr. Executive QA	Quality Assurance	To prepare process validation protocol and Report
4.	Production Manager	Production Manager	Overall responsible for checking & implementation of process validation
5.	QC Executive/ Manager	Quality Control	To analyze all the sample drawn at different stages during process validation
6.	QA Head	Quality Assurance	Overall responsible for compliance of Process Validation
7	Quality Head	Quality Assurance	Responsible for approval of Process Validation Protocol & report



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 6 Of 17

7. Product & Batch Detail:

Process Validation of **Lenalidomide Capsules 10 mg** shall be carried out on any three consecutive batches.

The details are as under:

Product/ Generic Name	:	Lenalidomide Capsules 10 mg
Strength	:	10 mg per capsules
Product Code	:	LD/003
Batch Size	:	500 Capsules
Dosage Form	:	Hard gelatin capsule for oral use only
Overage	:	2 %
Composition	:	Each hard gelatin capsule contains: Lenalidomide 10 mg Excipients q.s Approved colour used in empty capsule shell
Packing Detail	:	white to off white colour granular powder filled in size “2” hard gelatin capsule having yellow body & yellow cap. Such 30 capsule packed in 30 ml HDPE bottle
Storage Condition	:	Store between 15°C to 30°C. Protect from light & moisture.

Batch Detail for validation

3 consecutive batches are selected

Sr. No.	Batch No.	Batch Size	Batch Manufacturing Date	Mfg. Date	Exp. Date
1.	OC0079	100000 Capsules	03/04/2024	04/2024	03/2027
2.	OC0080	100000 Capsules	05/04/2024	04/2024	03/2027
3.	OC0081	100000 Capsules	07/04/2024	04/2024	03/2027



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 7 Of 17

8. Equipment and machinery List

Sr. No.	Name of equipment / instruments	ID No.
1.	Granulation Isolator	KPL/CT/001
2.	Automatic Capsule Filling Machine with line machine	KPL/CT/003
3.	Blister Packing	KPL/CT/004
4.	Octagonal Blender	KPL/CT/008
5.	Octagonal Blender	KPL/CT/009
6.	Starch Paste Kettle	KPL/CT/010
7.	Vibro shifter	KPL/CT/011
8.	Rapid Mixer Granulator	KPL/CT/012
9.	Fluid bed dryer	KPL/CT/013
10.	Cone mill	KPL/CT/014
11.	Colloid mill	KPL/CT/015
12.	Metal detector	KPL/CT/016
13.	Electronic Weighing Balance	KPL/CT/017
14.	Electronic Weighing Balance	KPL/CT/018
15.	Electronic Weighing Balance	KPL/CT/019
16.	Electronic Weighing Balance	KPL/CT/022
17.	Fogger	KPL/CT/023
18.	Leak Test	KPL/CT/024
19.	Electronic Weighing Balance	KPL/CT/026
20.	Electronic Weighing Balance	KPL/CT/027



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 8 Of 17

9. Batch Weighing Sheet:

Batch Size: 500 Capsules

Sr. No	Name of Ingredient	Spec.	Qty per capsule (in mg)	Overage (%)	Qty for batch 1.0 Lac (in Kg)
1	Lenalidomide	IHS	10.000	2%	1.0400 kg
2	Croscarmellose Sodium	BP	12.500	Nil	1.250 kg
3	Microcrystalline cellulose-102	BP	80.000	Nil	8.000 kg
4	Lactose	BP	173.550	Nil	17.3550 kg
5	PVK K-30	BP	2.850	Nil	0.285 kg
6	Isopropyl Alcohol IPA	BP	50.000	Nil	5.000 kg
7	Talcum	BP	2.850	Nil	0.285 kg
8	Magnesium Stearate	BP	2.850	Nil	0.285 kg
9	Empty hard gelatin Capsule shell size # 2 having yellow body and yellow Cap.	IHS	1NOS	Nil	102000 Nos.

Average Fill Content of Capsules: 285 mg \pm 10 %

Average Fill Content of 20 Capsules: 5.7 gm \pm 2 %

Average Weight of Capsule: 347 mg \pm 10 %
(Filled Weight + Empty Capsule)



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 9 Of 17

10. Approved Raw and Packing Material:

The raw or packing material uses in manufacturing process are from approved vendor.

Raw Material:

Raw Material	Spec.	Vendor
Lenalidomide	IHS	Adley Lab Limited
Croscarmellose Sodium	BP	Avon Pharma Chemicals
Microcrystalline cellulose- 102	BP	Himanshu Traders
Lactose	BP	Avon Pharma Chemicals
PVP K-30	BP	Anand Lathiya
Isopropyl Alcohol IPA	BP	Avon Pharma Chemicals
Talcum	BP	Neelkanth Minechem
Magnesium Stearate	BP	Zacfa Chemicals

Packing Material:

Packing Material	Spec.	Vendor
30 ml Flat HDPE Bottle	IHS	Preet Plastics, INDIA
Printed Sticker Label	IHS	Arihant Labels
Printed Laminated Carton	IHS	Sain Packaging Pvt. Ltd.




KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 10 Of 17

11. PROCESS FLOW DIAGRAM





KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 11 Of 17

12. Manufacturing Process

Check the Equipment for Cleanliness prior to use before manufacturing process and check the status label displayed. Take line clearance from IPQA for cleanliness

12.1 Sifting:

Sifting process carried out to obtain granules of uniform size and to remove extrinsic particles and this having impact on flow and filling and other physical properties of Capsules.

Sift the raw material as per the sieves size mentioned below

Sr. No	Raw Material	Sieve Size
1.	Lenalidomide	60 #
2.	Croscarmellose Sodium	60 #
3.	Microcrystalline cellulose-102	40 #
4.	Lactose	40 #
5.	PVK K-30	40 #
6.	Isopropyl Alcohol IPA	40 #
7.	Talcum	40 #
8.	Magnesium Stearate	40 #

Check the integrity and cleanliness of sieves before and after sifting process.

Sr. No	Sieve No	Check Sieve Integrity before Sifting	Checked by		Check Sieve Integrity after Sifting	Checked by	
			Production	QA		Production	QA
1.	40 #						
2.	60 #						

12.2 Dry Mixing:

1. Transfer Lenalidomide, Microcrystalline Cellulose into blender and mix for 10 min.
2. Transfer the Croscarmellose Sodium, Lactose and PVPK-30 & Talcum in step 1. And mix for specified time and visually examined the ingredients; it should be uniformly mixed and dense

Sr. No	Mixing Time (minutes)
1.	5 Minutes
2.	10 minutes
3.	15 minutes



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 12 Of 17

IPQA personal shall draw 1.0 gm of Samples each at different mixing time and send to Quality Control Department for Analysis

12.3 Lubrication:

Unload the previously mixed granules in Double cone blender and add magnesium Stearate. Mix the whole blend for specified time & visually examine the blend; the blend should be properly mixed and uniformly dense.

Sr. No	Mixing Time (minutes)
1.	3 minutes
2.	5 minutes

IPQA personal shall draw 1.0 gm of Samples each at different mixing time and send to Quality Control Department for Analysis.

Unload the blend sample in a double lined polythene bag placed in HDPE drum, close the polythene tightly and transfer the bulk to Capsule filling area.



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 13 Of 17

13. Filling & Polishing:

Take line clearance from IPQA for cleanliness of area and equipment. Record the temperature and humidity of the Capsule Filling Area

Assemble the capsule filling machine and put the blend in the hopper and run the capsule filling machine at different speed and check the physical parameter. Discard the capsules obtained during trial run.

Sr. No	Speed
1.	Minimum Speed
2.	Optimum Speed
3.	Maximum Speed

Inform IPQA personal for sampling of Filled capsules.

IPQA personnel shall draw sample of Filled Capsules at minimum speed, optimum speed & maximum speed and send to Quality Control Department for Analysis.

Load the filled capsule in the polishing machine to remove the debris of powder sticking with the capsule

14. Visual Inspection

Visual inspection helps to remove the physical defects of the capsules, thus ensure appearance quality of Capsules.

Inspect the filled and polished capsules for any defects.

Sr. No	Physical Defects/Rejections
1.	Empty Capsules
2.	Splitting
3.	Cracking
4.	Denting
5.	Mottling



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 14 Of 17

15. Sampling Plan:

S. No.	Stage	Sampling time / Qty.	Test / Process parameters
1	Mixing	After Mixing of 5 minute 10 minute 15 minute Sample quantity approx. 1.0 gm for each sample	Appearance Assay
2	Lubrication	After Mixing of 3 minute 5 minute Sample quantity approx. 1.0 gm for each sample	Appearance Assay
3	Capsule Filling & Polishing	Speed: Minimum Speed Optimum Speed Maximum Speed Sample quantity approx. 20 capsule	Description / Appearance Average Fill Weight of capsule Uniformity of weight/Content DT / Dissolution Assay
4.	Visual Inspection	At different speed of capsule filling machine	Empty Capsules Splitting Cracking Denting Mottling
5.	Packing	Cleanliness of bottle Extrinsic particles (Dust or black particles) Leak test	Should be clean and dry Should be free from extrinsic particles Should be comply



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 15 Of 17

16. Acceptance Criteria

Stage	Test / Process Parameter	Acceptance Criteria
Dry Mixing	Appearance	When observed visually material should appear uniform and dense.
	Assay	90.0 % to 110 %
Lubrication	Appearance	When observed visually material should appear uniform and dense
	Assay	90.0 % to 110 %
Filling & Polishing	Appearance	white to Off white colour granular powder filled in size "2" hard gelatin capsule having yellow body & yellow cap
	Average fill content	285 mg \pm 10 % (256.5 mg to 315.5 mg)
	Average weight of filled capsule	347 mg \pm 10 % (312.3 mg to 381.7 mg)
	Uniformity of dosage Units (By Weight Variation)	Not more than two of the individual weights deviate from the average weight by more than the 10 percent and none should deviate more than 20 %
	Dissolution	Not Less Than 75 % (Q) of Lenalidomide dissolved in 45 minutes
	Assay	90 % to 110 %
Visual Inspection	Empty Capsules	Should be Nil, or not more than 1.5 % of total batch
	Splitting	
	Cracking	
	Denting	
	Mottling	
Packing	Cleanliness of bottle	Should be clean and dry
	Extrinsic particles (Dust or black particles)	Should be free from extrinsic particles
	Leak test	Should be comply



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 16 Of 17

17. Deviation Report:

18. Change Control:

19. Abbreviation:

PVP	:	Process Validation Protocol
IPQA	:	In-process Quality Assurance
USP	:	United State Pharmacopoeia
BP	:	British Pharmacopoeia
IHS	:	In House Specification
QA	:	Quality Assurance
mg	:	milligram
gm	:	gram
q.s	:	quantity sufficient
mm	:	millimeter
AC	:	Anticancer
Hr	:	Hour
S.S	:	Stainless Steel
HDPE	:	High density polyethylene



KWALITY PHARMACEUTICALS LTD.
1-A, INDUSTRIAL AREA, RAJA KA BAGH, TEH. NURPUR, DISTT.
KANGRA (H.P) 176201 INDIA

PROCESS VALIDATION PROTOCOL

Product Name	Lenalidomide Capsule 10 Mg	
Protocol No.	LD/PVP/004-00	Supersedes No.: Nil
Effective Date	03/04/2024	Page 17 Of 17

20. Post Approval

Protocol Checked By

Name	Designation	Signature	Date

Protocol Approved By

Name	Designation	Signature	Date